



Gustafson 

Batch Lab Treater (BLT)

Recommended for Scale-up Applications and Laboratory Samples

FEATURES & BENEFITS

Designed to accurately apply seed treatment products in the laboratory to all seed types, for start-up treating runs. Engineered for quick treatment changes.

Caution: This treater is designed for laboratory lots!
DO NOT use for commercial or continuous operation.

- Light weight, portable and durable
- Internal lift bars help mix treatment with seed
- Bowl adjustment arm with tilt-lock angle
- Injection and compressed air provide accurate distribution of treatment-to-seed coverage
- Stainless steel mixing bowls: 1/4 & 5 pounds (120g & 2.2kg)

More Features & Benefits on back





Treater Dimensions: Height 28.69" Length 22.58" width 12.02"



Bowl Gear Motor

- 1/15hp, 1ph, 115V, TEFC, 50/60Hz
- 115V Power cord



Chemical Applicator Air Compressor

- 115V, 1ph, 50 or 60Hz (specify)
- 115V Power cord



Chemical Applicator

- Gun assembly with graduated, disposable syringe



Controls

- Bowl timer, 5 minute, self-contained and interlocked
- ON/OFF pump & bowl switches



Rubber Foot Mounts

- Absorbs vibration, scratch-free



Bowl Tilt Arm

- Adjustable (bowl angle)
- Locking mechanism

Value Added Features

- ✓ 2 mixing bowl sizes



- ✓ Mixing bowl lift bars



IMPORTANT: Technical information expressed in this document is for reference purposes only and subject to change without notice. If certified dimensional drawings are required call 1.800.634.6738